english

MiniSpacer® Dual Spray Metered Dose Inhaler Adapter

Bi-Directional Spray and Cap
22 mm O.D. / 22 mm I.D. ISO Connections

Accepts plastic and metal pMDI canister stems, including canisters with integrated dose counters.

INDICATIONS FOR USE

Model number 1025 is a single patient, disposable device for dispensing pressurized metered dose inhaler (pMDI) medication into a breathing circuit, as prescribed by a physician or other licensed health care practitioner. The device is indicated for patients on a breathing circuit, for whom aerosol medication has been prescribed, in short and long term critical care environments. REF 1025 is intended for use only when connected to a 22 mm fitting.

DIRECTIONS FOR USE

Installing the MINISPACER Adapter in the Patient Circuit

1. Remove the MINISPACER adapter from package.
2. Connect the MINISPACER device to the inspiratory limb of the patient circuit at the patient wye or according to facility guidelines.

CAUTION : Ensure the device is securely installed before use (i.e. the patient wye is fully inserted into the MINISPACER 22 mm I.D. connector and the patient tubing is securely over 22 mm O.D. connector).

NOTE : When correctly installed, the MINISPACER MDI port will be vertical and the pMDI canister will point downward during treatment.

3. Place the cap firmly on the MINISPACER port.

Administering MDI medication via the MINISPACER Adapter

1. Uncap the MINISPACER MDI port.
 2. Insert the MDI canister stem into the MINISPACER port.
 3. Coordinate MDI actuation with the beginning of inspiration.
- CAUTION :** Support the MINISPACER adapter while depressing the canister.
4. Administer the dose as ordered.
 5. After medication has been delivered as prescribed, remove the pMDI canister and firmly cap the MINISPACER port.
 6. Monitor the patient response.

NOTE : Discard the MINISPACER adapter when replacing the circuit according to facility guidelines, not to exceed two weeks.

CAUTION: REF 1025 MINISPACER adapter will not increment the counter on pMDI canisters with integrated dose counters surrounding the stem.

WARNING : Do not clean and reuse or sterilize. Potential risks associated with cleaning and reusing, reprocessing and/or sterilization include: cross-contamination or infection resulting from inadequate cleaning, product degradation, and/or failure of the device to perform as intended.

Symbol Glossary

	Manufacturer
	Use by Date
	LOT
	CE
	REF

	Consult instructions for use
	Authorized representative in the European Community
	Prescription only. U.S. Federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner.
	Do Not Reuse

Item No.	Description
REF 1025	MiniSpacer® Dual Spray MDI Adapter, 22 mm O.D. / I.D.

Thayer Medical Corporation
4575 S. Palo Verde Rd.
Tucson, AZ 85714-1961 USA
www.thayernedical.com
Made in USA

0482

G-1541 Rev. D

* Trademarks and registered trademarks of Thayer Medical Corporation.
© Thayer Medical Corporation 2015. All rights reserved.

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** trademarks and registered trademarks of their respective companies.

The MINISPACER dual-spray adapter was evaluated for particle size and dosing characteristics using an 8 stage cascade impactor. The MINISPACER was connected to the cascade impactor with ventilator tubing, wye and an endotracheal (ET) tube as they would be used in a ventilator circuit. Air was drawn through the assembly at a flow rate of 28.3 l/min. The aerosol samples exited the ET tube into the cascade impactor and comprised the Total Delivered Dose.

REF 1025 22 mm O.D. / I.D. – Installed 6 inches upstream of the wye in the inspiratory limb of a 22 mm circuit

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