MiniSpacer\* Dual Spray Metered Dose Inhaler Adapter With Counter Incrementing Actuator

Bi-Directional Spray and Cap

22 mm O.D. / 22 mm I.D. ISO Connections

**INDICATIONS FOR USE**

Model number 1025A is a single patient, disposable device for dispensing pressurized metered dose inhaler (pMDI) medication into a breathing circuit, as prescribed by a physician or other licensed health care practitioner. The device is indicated for patients on a breathing circuit, for whom aerosol medication has been prescribed, in short and long term critical care environments. REF 1025A is intended for use only when connected to a 22 mm fitting.

**DIRECTIONS FOR USE**

**Installing the MINISPACER Adapter in the Patient Circuit**

- Remove the MINISPACER adapter from package.
- Connect the MINISPACER device to the inspiratory limb of the patient circuit at the patient wye or according to facility guidelines. CAUTION : Ensure the device is securely installed before use (i.e. the patient wye is fully inserted into the MINISPACER 22 mm I.D. connector and the patient tubing is securely over 22 mm O.D. connector). NOTE : When correctly installed, the MINISPACER MDI port will be vertical and the pMDI canister will point downward during treatment.
- Place the cap firmly on the MINISPACER port.

**Administering MDI medication via the MINISPACER Adapter**

- Uncap the MINISPACER MDI port.
- Insert the MDI canister stem into the MINISPACER port.

NOTE : If canister stem features a dose-counting mechanism in the tip, align the MINISPACER's gear with the dose-counter. When correctly aligned, the canister's number dial will be visible through the MINISPACER window.

- Coordinate MDI actuation with the beginning of inspiration.

**CAUTION : Support the MINISPACER adapter while depressing the canister.**

- Administer the dose as ordered.
- After medication has been delivered as prescribed, remove the pMDI canister and firmly cap the MINISPACER port.
- Monitor the patient response.

NOTE : Discard the MINISPACER adapter when replacing the circuit according to facility guidelines, not to exceed two weeks.

WARNING : Do not clean and reuse or sterilize. Potential risks associated with cleaning and reusing, reprocessing and/or sterilization include: cross-contamination or infection resulting from inadequate cleaning, product degradation, and/or failure of the device to perform as intended.

### Symbol Glossary

	Manufacturer		Consult instructions for use
	Use by Date		Authorized representative in the European Community
	Batch code/lot number		Prescription only. U.S. Federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner.
	CE Mark		Do Not Reuse
	Catalog, reorder or reference number		Caution. The symbol can also be used to mean Attention, see Instructions for Use.

Item No.	Description
 1025A	MiniSpacer® Dual Spray MDI Adapter, 22 mm O.D. / I.D.

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G-1587 Rev. D

\* Trademarks and registered trademarks of Thayer Medical Corporation.

© Thayer Medical Corporation 2015. All rights reserved.

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<sup>[\*]</sup> trademarks and registered trademarks of their respective companies.

The MINISPACER dual-spray adapter was evaluated for particle size and dosing characteristics using an 8 stage cascade impactor. The MINISPACER was connected to the cascade impactor with ventilator tubing, wye and an endotracheal (ET) tube as they would be used in a ventilator circuit. Air was drawn through the assembly at a flow rate of 28.3 l/min. The aerosol samples exited the ET tube into the cascade impactor and comprised the Total Delivered Dose.

REF 1025A 22 mm O.D / I.D. Installed at the patient wye in the inspiratory limb of a 22 mm circuit.

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