1 LiteAire® Basic

THAYER MEDICAL CORPORATION  
4575 S. Palo Verde Rd., Ste 337  
Tucson, AZ 85714 USA  
www.ThayerMedical.com


CE

Rx only

MD

EC REP

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\*Below Limit of Quantification (BLOQ)    \*\*Trademarks and registered trademarks of their respective companies

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References/Referencias/Références/Referenzen/Bibliografia/Referenties:

1.

Lavorini et al. Spacers and Valved Holding Chambers—The Risk of Switching to Different Chambers. J. Allergy Clin. Immunol. Pract. 2020; 8(5): 1569-1573. Spacers and Valved Holding Chambers—The Risk of Switching to Different Chambers - ScienceDirect

2.

Mclvor et al. Optimizing the Delivery of Inhaled Medication for Respiratory Patients: The Role of Valved Holding Chambers. Can. Respir. J. 2018; 2018: 1-8. Optimizing the Delivery of Inhaled Medication for Respiratory Patients: The Role of Valved Holding Chambers (hindawi.com)

3.

Vincken et al. Spacer devices for inhaled therapy: why use them, and how? Eur. Respir. J. 2018; 4(2): 00065-2018. Spacer devices for inhaled therapy: why use them, and how? | European Respiratory Society (ersjournals.com)

4.

Johnson et al. The effect of a holding chamber on albuterol metered-dose inhaler product differences. Ann. Allergy Asthma Immunol. 2016; 117:246-250.

5.

Levy et al. Asthma patients' inability to use a pressurized metered-dose inhaler (pMDI) correctly correlates with poor asthma control as defined by the global initiative for asthma (GINA) strategy: a retrospective analysis. Prim Care Respir J 2013; 22:406-411. Asthma patients' inability to use a pressurised metered-dose inhaler (pMDI) correctly correlates with poor asthma control as defined by the Global Initiative for Asthma (GINA) strategy: a retrospective analysis | Primary Care Respiratory Journal (nature.com)

REF

1803 LITEAIRE® BASIC VALVED HOLDING CHAMBER

EN

## INSTRUCTIONS FOR USE (single-patient use only)

### INTENDED USE / INTENDED PURPOSE

The LiteAir® is a collapsible, disposable, dual-valved holding chamber designed to aid in the delivery of aerosolized medications delivered via a pressurized metered dose inhaler (MDI). The LiteAir® features a standard port designed for compatibility with standard MDI mouthpieces. It is a non-sterile device for single-patient use.

### INDICATIONS FOR USE

The LiteAir® is intended to be used by adults, adolescents and children ages 5 and up who are able to use a holding chamber without the aid of a mask and who are under the care or treatment of a physician or licensed healthcare professional. The environment for use includes homes, hospitals, and clinics.

### CONTRAINDICATIONS

The LiteAir® should not be prescribed for individuals who require the use of a mask.

### EXPECTED CLINICAL BENEFITS

The LiteAir® is a valved holding chamber, which is a device designed to reduce oropharyngeal deposition of large particle aerosolized drugs without compromising small particle drug delivery to the lungs<sup>1,2,3,4</sup>. Valved Holding Chambers also reduce the need to coordinate MDI actuation with inhalation<sup>1,2,3,4</sup>. Compared to MDIs used alone, the use of a Valved Holding Chamber can improve airway hyperresponsiveness<sup>1</sup>, lung function<sup>1</sup>, asthma control<sup>1,5</sup>, and reduce the local<sup>1,2,3</sup> and systemic<sup>1,4</sup> side effects of inhaled corticosteroids.

### HOW TO USE YOUR LITEAIRE®

The LiteAir® can be used right out of the package. Prior to use, ensure these instructions and the instructions supplied with the MDI have been read. Always follow your physician's instructions.

- Open packaging and remove LiteAir® device. (Fig. 1)
- Inspect the LiteAir® to ensure there is no damage to the device such as paperboard not adhered, loose valves, tears, rips, holes, discoloration, etc. Do not use the LiteAir® if device or packaging is damaged or compromised.
- Place your fingers in the middle of the LiteAir® along the long sides of the chamber. (Fig. 2)
- Gently squeeze the device to pop up the device and snap the side panels into place. (Fig. 2)
- Press inward where the inhaler port is located until it turns concave. (Fig. 3)
- Remove the cap from the inhaler mouthpiece. Consult the MDI manufacturer's instructions to shake and prime the inhaler prior to use (Fig. 4)
- Look inside the mouthpiece and chamber of the LiteAir®. Remove any small objects that could be inhaled. (Fig. 5)
- Insert the inhaler mouthpiece firmly into the LiteAir® inhaler port. (Fig. 6)
- Firmly squeeze your LiteAir®, pressing the fingerprint icons at the point where you feel the cardboard inner wall, to ensure the most effective treatment. Close your lips firmly around the small holes on the mouthpiece of the LiteAir® and exhale fully. (Fig. 7)
- If you inhale, spray one puff of medication into the chamber. Continue inhaling slowly and deeply. Hold your breath for 5 to 10 seconds, or as directed by the MDI manufacturer or your physician. (Fig. 8)
- If your physician has directed you to take more than one puff, repeat steps 9-10.
- Remove the inhaler from the inhaler port. (Fig. 9)
- Collapse the LiteAir® by pressing gently until flat. (Fig. 10)

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### STORAGE AND HANDLING INSTRUCTIONS

- Store your LiteAir® in a clean, dry location inside an unsealed bag (Fig. 10). Do not store your LiteAir® in a sealed or airtight container as this will not allow the LiteAir® to dry between uses.
- Each LiteAir® is good for multiple uses by the same user. Do not share your LiteAir® with anyone.
- Replace your LiteAir® at least weekly, or sooner if it becomes damaged, crushed, soggy, contaminated or discolored. The LiteAir® may be placed in the general household waste receptacle.
- If LiteAir® devices from multiple patients are stored together, ensure they do not come into direct contact with each other to prevent the potential exchange of pathogens. Note: If sharing an MDI between parties (shared canister protocol) ensure the MDI is thoroughly cleaned/disinfected between uses.

### ⚠ LIMITATIONS, CAUTIONS, PRECAUTIONS OR WARNINGS

#### WARNING

Always look inside the mouthpiece and chamber before you use the LiteAir® to be sure there are no small objects that could be inhaled (such as coins, lint, buttons, etc.)

#### CAUTION

Do not clean, disinfect or sterilize your LiteAir®. The LiteAir® is made from paper products and may be permanently damaged if cleaned, disinfected or sterilized. Potential risks include: product degradation and failure of the device to perform as intended.

#### CAUTION

The device is not anticipated to cause any undesirable side effects; however, if the user experiences undesirable side effects from the administered MDI medication, they should consult with a physician immediately. Please refer to the MDI medication Instructions for Use for more information. If symptoms persist, seek immediate medical treatment.

#### CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician or other practitioner licensed to prescribe by law.

#### CAUTION

The LiteAir® should not be used with the following: dry powder inhalers, breath-actuated inhalers, or soft mist inhalers. Refer to the instructions supplied with your inhaler for more information.

#### NOTICE OF SERIOUS INCIDENTS

Please report any suspected serious incidents in relation to this device to the manufacturer using the contact information below, and contact the competent authority of the Member State in which the incident occurred.

#### CONTACT INFORMATION

To obtain information regarding LiteAir® and its uses or to report a suspected serious incident, please contact your local representative or Thayer Medical directly at 1-800-250-3330, email at [info@thayermedical.com](mailto:info@thayermedical.com), or visit our website at [www.thayermedical.com](http://www.thayermedical.com).

### THE LITEAIRE® BASIC VALVED HOLDING CHAMBER PARTICLE SIZE AND DOSING CHARACTERISTICS

The LiteAir® Valved Holding Chamber was evaluated for particle size and dosing characteristics using an 8-stage cascade impactor. The LiteAir® was connected to the induction port of the cascade impactor with a mouthpiece adapter. Air was drawn through the assembly at a flow rate of 28 L/min (simulating adult use) and 14 L/min (simulating pediatric use). Three LiteAir® were used 3 times each to evaluate each drug (N=9 evaluations/drug). See table on reverse for results.

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