MiniSpacer* Dual Spray Metered Dose Inhaler Adapter With Counter Incrementing Actuator

Bi-Directional Spray and Cap

22 mm O.D. ISO Connections

INDICATIONS FOR USE

Model number 1024A is a single patient, disposable device for dispensing pressurized metered dose inhaler (pMDI) medication into a breathing circuit, as prescribed by a physician or other licensed health care practitioner. The device is indicated for patients on a breathing circuit, for whom aerosol medication has been prescribed, in short and long term critical care environments. REF 1024A is intended for use only when connected to a 22 mm fitting.

DIRECTIONS FOR USE

Installing the MINISPACER Adapter in the Patient Circuit

- Remove the MINISPACER adapter from package.
- Connect the MINISPACER device to the inspiratory limb of the patient circuit, approximately 6 inches back from the patient wye or according to facility guidelines.

CAUTION : Ensure the tubing is secure over MINISPACER connectors when in use.

NOTE : When correctly installed, the MINISPACER MDI port will be vertical and the pMDI canister will point downward during treatment.

- Place the cap firmly on the MINISPACER port.

Administering MDI medication via the MINISPACER Adapter

- Uncap the MINISPACER MDI port.
- Insert the MDI canister stem into the MINISPACER port.

NOTE : If canister stem features a dose-counting mechanism in the tip, align the MINISPACER’s gear with the dose-counter. When correctly aligned, the canister’s number dial will be visible through the MINISPACER window.

- Coordinate MDI actuation with the beginning of inspiration.


CAUTION : Support the MINISPACER adapter while depressing the canister.

- Administer the dose as ordered.
- After medication has been delivered as prescribed, remove the pMDI canister and firmly cap the MINISPACER port.
- Monitor the patient response.

NOTE : Discard the MINISPACER adapter when replacing the circuit according to facility guidelines, not to exceed two weeks.

WARNING : Do not clean and reuse or sterilize. Potential risks associated with cleaning and reusing, reprocessing and/or sterilization include: cross-contamination or infection resulting from inadequate cleaning, product degradation, and/or failure of the device to perform as intended.

Symbol Glossary

| | | | |
|------------------------------------------------------------------------------------|--------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
|  | Manufacturer |  | Consult instructions for use |
|  | Use by Date |  | Authorized representative in the European Community |
|  | Batch code/lot number |  | Prescription only. U.S. Federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner. |
|  | CE Mark |  | Do Not Reuse |
|  | Catalog, reorder or reference number |  | Caution. The symbol can also be used to mean Attention, see Instructions for Use. |

| Item No. | Description |
|------------------------------------------------------------------------------------------|-------------------------------------------------------|
|  1024A | MiniSpacer® Dual Spray MDI Adapter, 22 mm O.D. / O.D. |

| | | | |
|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------|
|  | Thayer Medical Corporation |  | 0482 |
|  | 4575 S. Palo Verde Rd. Tucson, AZ 85714-1961 USA www.thayermedical.com Made in USA |  | MDSS GmbH Schiffgraben 41 30175 Hannover, Germany |

G-1586 Rev. D

* Trademarks and registered trademarks of Thayer Medical Corporation.

© Thayer Medical Corporation 2015. All rights reserved.

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^[**] trademarks and registered trademarks of their respective companies.

The MINISPACER dual-spray adapter was evaluated for particle size and dosing characteristics using an 8 stage cascade impactor. The MINISPACER was connected to the cascade impactor with ventilator tubing, wye and an endotracheal (ET) tube as they would be used in a ventilator circuit. Air was drawn through the assembly at a flow rate of 28.3 l/min. The aerosol samples exited the ET tube into the cascade impactor and comprised the Total Delivered Dose.

REF 1024A 22 mm O.D. / O.D. – Installed 6 inches upstream of the wye in the inspiratory limb of a 22 mm circuit

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