




Collapsible MDI Holding Chamber: An Innovative Alternative

LiteAire®'s unique dual-valved MDI holding chamber design delivers pop-up convenience and effective drug output at a fraction of the cost.

In most clinical settings, the LiteAire MDI holding chamber can reduce costs by replacing existing rigid plastic holding chambers or inefficient spacers with a paperboard alternative. This unique paperboard design allows the LiteAire to be reused by a patient over multiple doses while meeting, and often exceeding, the performance of plastic holding chambers.

LiteAire cross-section with the side panel removed.



Convenient, individually packaged, 25-count dispenser box.

- Dual Valved
- Pop Up To Use
- Portable/Stores Flat
- Labeled For One Week's Use
- Not Made With Natural Latex
- Single Patient Use

The LiteAire delivers true value when used in:

- **Emergency Room:** For patients requiring the use of a holding chamber for immediate treatment; ideal for isolation room conditions.
- **Inpatient/Outpatient Pulmonary Function Testing Lab:** Clean, disposable, ready to use, designed to accept most pressurized metered dose inhalers.
- **Short-Term Respiratory Infections:** A cost effective treatment method that meets the short-term needs of the acute respiratory patient.

INFORMATION / SAMPLES

800-250-3330
info@thayermedical.com

ORDERING INFORMATION

e-mail: orders@thayermedical.com
Fax: 520-790-5854

Single 25-Count Dispenser Box - Part No. 1304 ... \$73.75
Case of four 25-Count Dispenser Boxes - Part No. 1306 ...\$295.00

* Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.



4575 S. Palo Verde Rd, Ste 337
Tucson, Arizona 85714-1961
www.thayermedical.com

