




Quake®: Mucus Clearance Device

The Quake® is an innovative and unique device that provides mucus clearance therapy in a portable, easy-to-use, handheld device.

The Quake®'s soft touch handle allows the patient to regulate rotation frequency and vary percussive pressure differentials within the lungs during both inhalation and exhalation. The Quake® provides optimum performance in any patient position and offers effective results due to its "Flow Independent" design. The Quake® truly allows for individualized therapy that meets each patient's mucus clearance needs (www.thayermedical.com/quake.htm).

Disassembled
Quake® view.



- Position-Independent
- Flow-Independent
- Patient Controlled by Rotation of Handle
- Percussion on Inhalation and Exhalation
- Portable
- Hand-Held
- Fits Standard 22mm Face Mask

Easy to clean: The Quake®'s simple, two-part design enables the user to easily and quickly disassemble and reassemble the Quake for cleaning. It can be cleaned with soap and water, boiled in water, or washed in a dishwasher. The Quake® can also be put in a Quick-Clean™ Micro Steam™ bag and steam cleaned in the microwave.*

The Quake® is designed to meet the demanding mucus clearance needs of a variety of different types of patients, including:

- Long-term, sputum-producing respiratory diseases
- Short-term, sputum-producing respiratory infections

Feel the power. Experience the sensation.

Ordering Information

Single Count (Part No. 1297)
Case of 5 (Part No. 1293)

For more information, call Thayer Medical at 800-250-3330 or e-mail us at sales@thayermedical.com



Turn Handle
While Inhaling
and Exhaling.



Thayer Medical
4575 S. Palo Verde Rd., Ste. 337
Tucson, Arizona 85714-1961
(800) 250-3330
www.thayermedical.com

©Copyright 2007, Thayer Medical Corporation. All rights reserved.
Quake® is a registered trademark of Thayer Medical Corporation. QK07-2 rev. A

* "Respiratory Equipment Cleaning and Disinfection Guidelines." The Children's Hospital-Denver, Headley, Aree, Dolan, Susan A., Lucas, Trent. 3 September 2004. November 2004 <http://www.thechildrenshospital.org/share/policy/HospitalChart.pdf>

Copyright 2004, The Children's Hospital Denver

Quick Clean™ and Micro-Steam™ are trademarks of Medela
Medela is a registered trademark of Medela