




Quake®: Airway Clearance Device

The Quake® is an innovative and unique device that provides mucus clearance therapy in a portable, easy-to-use, handheld device.

The Quake®'s soft touch handle allows the patient to regulate rotation frequency and vary percussive pressure differentials within the lungs during both inhalation and exhalation. The Quake® provides optimum performance in any patient position and offers effective results due to its "Flow Independent" design. The Quake® truly allows for individualized therapy that meets each patient's mucus clearance needs

- Position-Independent
- Flow-Independent
- Patient Controlled by Rotation of Handle
- Percussion on Inhalation and Exhalation
- Portable
- Hand-Held
- Fits Most 22mm Connections

Disassembled Quake® view.



Easy to clean: The Quake®'s simple, two-part design enables the user to easily and quickly disassemble and reassemble the Quake for cleaning. It can be cleaned with soap and water and/or boiled in water.

The Quake® is designed to meet the demanding mucus clearance needs of a variety of different types of patients, including:

- Long-term, mucus-producing respiratory diseases
- Short-term, mucus-producing respiratory infections

Feel the power. Experience the sensation.

Ordering Information

Part Number 1297 - Available in cases of 5 eaches

For more information, call Thayer Medical at 800-250-3330 or e-mail us at sales@thayermedical.com




 Thayer Medical
 4575 S. Palo Verde Rd., Ste. 337
 Tucson, Arizona 85714-1961
 (800) 250-3330
www.thayermedical.com

© Copyright 2016, Thayer Medical Corporation. All rights reserved.

Quake® is a registered trademark of Thayer Medical Corporation. QK07-6 rev. B